



Handheld Ultrasound
Scanner Manufacturer



WIRELESS PROBE ULTRASOUND

AFFORDABLE
VISIBLE - FAST - EASY

GUANGZHOU SONOHEALTH MEDICAL TECHNOLOGIES CO., LTD.

Address: Room 601, Building 3, No.7 Ruitai Road, Huangpu District, Guangzhou, China

E-mail: hello@sono-health.com

Telephone: +86-02-87831007

Phone/WhatsApp: +86-15013156300

Website: www.sono-health.com



YOU CAN USE IT ANYTIME & ANYWHERE



APPLICATION

Application Sites: Abdomen, OB/GYN, Lung, Urology

Application Directions

- Emergency inspection: ER, ICU, Wild First Aid, Battle Field Rescue.
- Preliminary examination: Ward Inspection, primary clinic examination, medical examination, health screening, home care, family planning, etc.
- Visualization tools: Invasive intervention guide, surgical and therapy guidance.
- Remote diagnosis, consultation, training: Works on smart phone or tablet, easy.



FEATURES

- Exceptional image quality, meet the needs of clinical ultrasound.
- Wireless connected, can be conveniently used in surgery without fixing of cables.
- And by using the disposable protecting cover, it can easily solve the sterilization issue of the probe.



SPECIFICATION

- Scan Angle/Length: 60
- Frequency: 3.5/5.0 MHz
- Scan Mode: B, B/M, Color, PW
- Display Depth: 90~305mm
- Power: Built-in lithium battery, wireless or USB charge 4h, work 5h
- Working System: iOS, Android
- WiFi Type: 802.11n / 2.4G / 5G
- Dimension: 150x60x20mm
- Weight: 200g



WIRELESS & USB CHARGE

- Both wireless charge and USB charge you can choose.
- Easy and connect with a wireless power bank or USB charger.



EASY TO USE

Simple operation, short learning curve for people without ultrasound experience.



WiFi CONNECT

- WiFi connect with your Android phone, iPhone or iPad.

HIGH QUALITY IMAGES

See The Difference Great Image Quality Makes

