



Handheld Ultrasound
Scanner Manufacturer



WIRELESS PROBE ULTRASOUND

AFFORDABLE
VISIBLE - FAST - EASY

GUANGZHOU SONOHEALTH MEDICAL TECHNOLOGIES CO., LTD.

Address: Room 601, Building 3, No.7 Ruitai Road, Huangpu District, Guangzhou, China

E-mail: hello@sono-health.com

Telephone: +86-02-87831007

Phone/WhatsApp: +86-15013156300

Website: www.sono-health.com



YOU CAN USE IT ANYTIME & ANYWHERE



APPLICATION

Application Sites: Small Parts, MSK, Nerve, Vascular access

Application Directions:

- Emergency inspection: ER, ICU, Wild First Aid, battle field rescue.
- Preliminary examination: Ward Inspection, primary clinic examination, medical examination, health screening, home care, etc.
- Visualization tools: Invasive intervention guide, surgical and therapy guidance.
- Remote diagnosis, consultation, training: Works on smart phone or tablet, easy.



FEATURES

- Exceptional image quality, meet the needs of clinical ultrasound.
- Wireless connected, can be conveniently used in surgery without fixing of cables.
- And by using the disposable protecting cover, it can easily solve the sterilization issue of the probe.



SPECIFICATION

- Scan Angle/Length: 40mm
- Frequency: 7.5/10 MHz
- Scan Mode: B, B/M, Color, PW
- Display Depth: 20~100mm
- Power: Built-in lithium battery, wireless or USB charge 4h, work 5h
- Working System: iOS, Android
- WiFi Type: 802.11n / 2.4G/5G
- Dimension: 150x50x20mm
- Weight: 200g



WIRELESS & USB CHARGE

- Both wireless charge and USB charge you can choose.
- Easy and connect with a wireless power bank or USB charger.



EASY TO USE

Simple operation, short learning curve for people without ultrasound experience.



WiFi CONNECT

WiFi connect with your Android phone, iPhone or iPad.

HIGH QUALITY IMAGES

See The Difference Great Image Quality Makes

